




## BCI Meeting - June 7<sup>th</sup> 2023



HappyFeat: an interactive and efficient BCI framework for clinical applications

Arthur Desbois,  
Inria Paris, ARAMIS team, Paris Brain Institute (ICM)


# What is “BCI”?




# What is “BCI”?




# What is “BCI”?




## Motor execution VS Rest



## Motor imagery VS Rest



## Spectral Power Decrease



Desynchronization effect  
(Pfurtscheller et al, 1999)

- **Features of interest (FOI)**


- Selecting adequate FOIs is a **crucial step** for BCI performance.
- After EEG signals acquisition, an **analysis phase** is needed to select best FOIs.

- ➔ Scientific softwares (i.e., MATLAB)
- ➔ Manual step, expertise needed

- If this analysis phase is too long, a lot can change in the meantime:

- EEG sensors impedances
- Subject's brain behaviors
- Subject's attention & motivation


- ➔ Signal characteristics might be very different between Acq/Training phase and Online phase...



The “analysis phase” involves many manipulations. An example:

- Setting up “feature extraction” scenarios in OpenViBE...
- Finding FOIs through visualization...
- Setting up & running training scenarios in OpenViBE...
- ... and maybe going again through those steps multiple times until “correct” features have been found, or to account for inter-run variability

➔ Tedious, error-prone,  
hard to achieve in a limited time





# HappyFeat - Main Concepts

## Python-based framework for facilitating ML pipelines

### Main focus:

making Feature Selection  
+ Classif. Training phases  
**easy & fast**

**Analyze your data,  
select your Features  
& train your classifier  
in less than 5 minutes!**

The screenshot displays the 'HappyFeat - Feature Selection interface' with three main sections:

- == FEATURE EXTRACTION ==**
  - Test files: Test-[2022.05.05-15.36.18].ov, Test-[2022.05.05-15.46.20].ov, Test-[2022.05.05-15.55.59].ov, Test-[2022.05.05-16.19.27].ov, Test-[2022.05.05-16.30.16].ov
  - Extraction parameters:
    - Epoch of Interest (EOI) (s): 3
    - EOI offset (s): 1
    - Sliding Window (Burg) (s): 0.25
    - Window Shift (Burg) (s): 0.161
    - Connectivity Metric: Magnitude Squared Coh.
    - Length of a connectivity measure (s): 0.25
    - Overlap btw. connectivity measures (%): 36
    - Auto-regressive estim. length (s): 0.038
    - Frequency resolution (ratio): 1
  - Experiment parameters (set in Generator GUI):
    - Nb Trials per class: 20
    - Class / Stimulation 1: MI
    - Class / Stimulation 2: REST
    - "Get Set" time (s): 20
    - Pre-Stimulus time (s): 3
    - Trial duration (s): 3
    - Inter-trial interval min (s): 2.5
    - Inter-trial interval max (s): 3.5
    - Feedback time (s) (online scenario): 3
  - Buttons:  Browse for OpenVIBE designer script, Extract Features and Trials
- == VISUALIZE FEATURES ==**
  - Test-[2022.05.05-16.30.16]-SPECTRUM-CONNECT
  - Figure 2: Topography of power spectra, for freq. 10:13 Hz. The figure shows a circular brain map with a color scale from 0.0 to 0.7.
  - Load file(s) for analysis:
    - Frequency min: 0
    - Frequency max: 40
    - Sensor for PSD visualization: CP3
    - Topography Freq (Hz), use ":" for freq band: 10:13
    - Scale Colormap (R²map and Topography): ☒
    - Buttons: Power Spectrum, Node Strength, Freq.-chan. R² map, Freq.-chan. R² map, PSD for the 2 classes, NodeStr. for the 2 classes, Brain Topography, Brain Topography
- == CLASSIFIER TRAINING ==**
  - Ex: FCz;14, FCz;14:22 (for freq range)
  - Power Spectrum: Add feature, Remove feature, CP3;10
  - Connectivity: Add feature, Remove feature, CP3;12
  - Number of k-fold for classification: 10
  - Test-[2022.05.05-16.30.16]-TRIALS.csv
  - Last Training Results:
    - Test-[2022.05.05-16.30.16]-TRIALS.csv
    - Feature(s) for PowSpectrum: Channel CP3 at 10 Hz
    - Feature(s) for Connectivity: Channel CP3 at 12 Hz
    - Overall accuracy : 100.0%
    - Class 1 | Precision : 1.0 | Sensitivity : 1.0 | F\_1 Score : 1.0
    - Class 2 | Precision : 1.0 | Sensitivity : 1.0 | F\_1 Score : 1.0
  - Buttons: TRAIN CLASSIFIER, FIND BEST COMBINATION





- **Key feats & mechanisms:**
  - **Clean, risk-free environment**
    - ➔ avoid unnecessary & error-prone manipulations.
  - **Trial-and-error oriented workflow**
    - ➔ all steps can be repeated quickly & as many times as needed
  - **Unified “dashboard” GUI**
  - **OpenViBE** used in the background, as a fast & efficient processing engine.
    - ➔ no scenario edition/manipulation necessary: everything is **automated!**
- **Two main use cases:**
  - **Make BCI pipelines smoother/easier to use** and allow **reproducibility** of exps.
  - **Prospective works & comparison of alternative features of interest** (connectivity, networks...)


- **Efficient processing pipelines**
  - Available features for classification: *Spectral Power*, *Connectivity-based network metrics*
  - ... It's also possible to train the classifier using **a fusion of both features**.
- **Feature extraction**
  - Easy access to all experiment & signal processing parameters.
  - Use **pre-recorded signals** or **on-the-fly** during acquisition phase.
- **Visual Analysis for feature selection**
  - $R^2$  maps, PSD comparison across trials, time/freq. ERD/ERS analysis, brain topography...
- **Classifier training**
  - **Run as many training attempts as needed**, using different features, in only a few clicks.
  - Concatenate trials from multiple recorded sessions
  - OpenViBE scenarios are **updated and launched on-the-fly**.
  - Automatically generates/updates the **“online classification”** scenario.

## 1. MI Pipeline / “Feature type” selection


- Selecting btw. multiple “template” scenarios depending on the use case (power spectrum, connectivity type...)
- Edit basic/common parameters (acquisition, extraction, training...)



Extraction scenario



Training scenario




Online classif.  
scenario


## 1. MI Pipeline / “Feature type” selection

## 2. Feature Extraction

- **Select signal files, and extraction parameters** (lengths and overlap of windows, FFT size...)
- **Run the generated extraction OpenViBE scenario (in the background)** for all selected signal files:
  - ➔ **Extract metadata** (sampling freq, electrodes...)
  - ➔ **Cut the signal to regions of interest** (MI trials & baseline portions), **generate CSV file** with only these chunks (for the training step)
  - ➔ **Apply a signal processing pipeline** (PSD computation, connectivity measure...) to the signal chunks of interest, **generate CSV files for future analysis**
- **Runs in an autonomous thread:**  
You can do visualizations and training attempts for signals already processed in the meantime.



Extraction scenario



Signal 1  
(.ov)

- ➔ Spectral Power (class 1) CSV
- ➔ Spectral Power (class 2) CSV
- ➔ Trials CSV
- ➔ Metadata CSV

Signal 2  
(.ov)

- ➔ Spectral Power (class 1) CSV
- ➔ Spectral Power (class 2) CSV
- ➔ Trials CSV
- ➔ Metadata CSV


...

1. MI Pipeline / “Feature type” selection
2. Feature Extraction
3. Analysis, Feature Selection

- Select one or multiple signals & load their spectral/connectivity data (CSV work files generated during “Feature Extraction”)
- Use different **Visualization Tools** to help find & select **features of interest (FOIs)** for training
  - ➔ Frequency/channel  $R^2$  map
  - ➔ PSD (or connect. metric) comparison btw. the 2 conditions (MI/REST) for a given electrode
  - ➔ Time/frequency ERD/ERS analysis for each condition
  - ➔  $R^2$  mapped as a brain topography for a given frequency (or range)




Combine as many visualization windows as you need

A “Dual metric” pipeline allows to show (for ex.)  $R^2$  maps for both Power Spectrum and Connectivity in parallel



Signal 1  
(.ov)

- ➔ Spectral Power (class 1) CSV
- ➔ Spectral Power (class 2) CSV
- ➔ Trials CSV
- ➔ Metadata CSV



1. ML Pipeline / “Feature type” selection
2. Feature Extraction
3. Analysis, Feature Selection
4. Classifier Training


- Set one or more **Features of Interest** (Electrode/Freq.)
- **Select file(s)** with which you want to train your classifier  
If > 1 file : their trials are automatically concatenated
- **Run the Classifier Training scenario**  
(Auto. generated in step 1 + auto. edited with FOIs)  
→ Classification **ACCURACY** + **WEIGHTS**

**Disappointing results?** (“My accuracy is 50% ☹”)  
Maybe try again with other features. It only takes a few seconds...

**Satisfying results?!** (“OMG 95%”)

Good news! The “Online Classification” scenario has already been automatically been updated with:

- Classifier training weights
- Features of interest used for training



Full list of dependencies:

- **Python 3.9**
  - shutils>=0.1.0
  - mne>=0.23.0
  - numpy>=1.21.1
  - pandas>=1.3.1
  - PyQt5>=5.15.6
  - statsmodels>=0.13.1
  - scipy>=1.7.1
  - spectrum>=0.8.0
  - matplotlib>=3.4.2
- **OpenViBE v3.5.0**



- **Current limitations**

- Low flexibility regarding electrode schemes. *(work in progress...)*
- Only one type of classification algo. proposed (LDA) *(work in progress...)*
- Pipelines are “fixed”:
  - trading OpenViBE’s high level of flexibility...
  - + ...for a high comfort of pipeline settings customization
  - + & “trial-and-error” workflow
- Only three types of pipelines/feature types:
  - Power Spectrum Density
  - Connectivity (coherence & its variants)
  - Dual (mixing PSD & Connectivity)
- ➔ In project (*long term!*) for more “**prospective power**”:
  - allow the user to choose btw. 1 and 3 feature-types & network metrics to mix as they see fit (MSC/node strength + Imag(Coh)/Laterality + ...)

- **Current limitations**

- Every time a new feature type is selected, or extraction parameters are modified...  
all work files need recomputating from scratch.
  - ➔ Necessity for a robust “work session” save/load mechanism
- In project: fully autonomous 100% Python version, without OpenViBE
  - No acquisition/online possibilities
  - Obviously slower... (no C++ optimizations!)
  - + More portable, all types of platforms supported (MacOS!)
  - + Other (more flexible) formats for I/O and work files (CSV, EDF...)

- **Already available online, work-in-progress version:**  
<https://github.com/Inria-NERV/happyFeat>
- **To be continued...**
  - **More flexibility** (pipeline options, mixing metrics...)
  - **More network metrics** (based on connectivity)
    - ... and associated visualization tools
  - **More options for classification algorithms**
  - **Workspace/session manager** to save/load session settings
  - **Fully autonomous Python version**, for cross-platform usage

**Stay tuned!**



*BCI Motor Imagery with OpenViBE in X-Men: First Class*

# Thanks for your attention! Any questions?